




Crystal Application Tutorial: Initial calculations


Load data




- Load the Crystal application (Click on the crystal application  on the top toolbar)
- Open the folder /RepTate/RepTate/data/Crystal/ and select the 3 files InitialTest***.shearxs



- Drag and drop all files on to the Reptate window





- You should now see




Run the model

- Load GO-polyStrand model (choose the model from the drop down list and then click


the theory icon ). The calculation should take about 4 seconds.



- The default parameters are fine for this calculation so click compute



- You should see




- Set the parameter N_0 $1e-5$ to produce agreement at very early time

<input checked="" type="checkbox"/>	tau0	3.8e-10	inf
<input checked="" type="checkbox"/>	rhoK	2.7e+09	inf
<input checked="" type="checkbox"/>	G_C	0.063	inf
<input checked="" type="checkbox"/>	N_0	1e-5	inf
<input type="checkbox"/>	beta	1	inf
<input type="checkbox"/>	delta	-0.5	inf

Things to try


- Change Gamma to adjust the sensitivity to shear
- Make *small* changes to epsilonB to adjust the quiescent barrier (equivalent to changing temperature)
- Adjust tau0 to scale all nucleation rates
- Adjust G_C to change the crystal growth rate

- Change the molecular weight distribution by clicking



and changing the relaxation times, number of modes, and concentrations (this needs to add up to 1).

- Double click on a data set and change the shear rate ($\dot{\gamma}$) and shear end time (t_{stop}) [T is temperature, which currently has no effect]



- Click the tabs on the left of the screen to make individual plots fill the window



- Use the 'View' menu towards the upper right corner to choose different views.

